





PDHA1 Antibody

Product Code	CSB-PA247984
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08559
Immunogen	Synthesized peptide derived from internal of Human PDHA1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PDHA1 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Pyruvate dehydrogenase E1 component subunit alpha; somatic form; mitochondrial; PDHE1-A type I; PDHA1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PDHA1

Image



Western blot analysis of extracts from A549, HepG2 and Jurkat cells, using PDHA1 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using PDHA1 antibody.