




RABEP2 Antibody

Product Code	CSB-PA953159
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H5N1
Immunogen	Synthesized peptide derived from internal of Human RABEP2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total RABEP2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	rab GTPase-binding effector protein 2; rabaptin-5beta; rabe2; rabpt5b;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RABEP2

Image



Western blot analysis of extracts from HeLa cells, Jurkat cells, 293 cells and HUVEC cells, using RABEP2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.