





RAB3GAP2 Antibody

Product Code	CSB-PA215524
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2M9
Immunogen	Synthesized peptide derived from internal of Human.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total RAB3GAP2 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	RAB3 GTPase activating protein 2 (non-catalytic); RAB3-GAP150; RAB3GAP150; RBGPR; RGAP-iso
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAB3GAP2

Image



Western blot analysis of extracts from COS cells, using RAB3GAP2 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using RAB3GAP2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.