




RABL2A Antibody

Product Code	CSB-PA967554
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBK7
Immunogen	Synthesized peptide derived from internal of Human RABL2A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total RABL2A protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Rab-like protein 2A;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RABL2A

Image



Western blot analysis of extracts from HuvEc and COLO205 cells, using RABL2A antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.