




RDX Antibody

Product Code	CSB-PA101431
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35241
Immunogen	Synthesized peptide derived from internal of Human RDX.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total RDX protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	RADI; RDX;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RDX

Image



Western blot analysis of extracts from HuvEc cells, using RDX antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.