




RANBP6 Antibody

Product Code	CSB-PA593893
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60518
Immunogen	Synthesized peptide derived from internal of Human RANBP6.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total RANBP6 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	RAN binding 6; RNBP6;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RANBP6

Image



Western blot analysis of extracts from K562 cells, COLO cells and COS-7 cells, using RANBP6 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.