




RAP1GAP Antibody

Product Code	CSB-PA918584
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47736
Immunogen	Synthesized peptide derived from C-terminal of Human RAP1GDS1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total RAP1GDS1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	KIAA0474; RAP1 GTPase activating protein; Rap1 GTPase-activating protein 1; RAP1; GTPase activating protein 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAP1GAP

Image



Western blot analysis of extracts from COLO205 cells, using RAP1GDS1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.