




RRAGA/RRAGB Antibody

Product Code	CSB-PA590465
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L523
Immunogen	Synthesized peptide derived from internal of Human RRAGA/B.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total RRAGA/B protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ras-related GTP-binding protein A; Rag A; RagA; Adenovirus E3 14.7 kDa-interacting protein 1; FIP-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RRAGA/RRAGB

Image



Western blot analysis of extracts from HepG2 cells, using RRAGA/B antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.