




RAB31 Antibody

Product Code	CSB-PA139202
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13636
Immunogen	Synthesized peptide derived from internal of Human RAB31.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total RAB31 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Ras-related protein Rab-31; Ras-related protein Rab-22B;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAB31

Image



Western blot analysis of extracts from HeLa cells and HepG2 cells, using RAB31 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.