



RAB34 Antibody

Product Code	CSB-PA087223
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZG1
Immunogen	Synthesized peptide derived from internal of Human RAB34.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total RAB34 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Rab-39; Ras-related protein Rah??ab-39; RAB39; RAH; ras-related protein Rab-34
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAB34

Image



Immunohistochemistry analysis of paraffin-embedded human skeletal muscle tissue using RAB34 antibody.