




RAB6C Antibody

Product Code	CSB-PA695641
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0N0
Immunogen	Synthesized peptide derived from internal of Human RAB6C.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total RAB6C protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	gtp-binding protein rab6c; rab gtpase rab6a; rab6; rab6a; member ras oncogene family
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAB6C

Image



Western blot analysis of extracts from COLO cells, using RAB6C antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.