





RGS1 Antibody

Product Code	CSB-PA274857
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08116
Immunogen	Synthesized peptide derived from internal of Human RGS1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total RGS1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	RGS1; Early response protein 1R20; B-cell activation protein BL34;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RGS1

Image



Western blot analysis of extracts from mouse muscle cells, using RGS1 antibody.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using RGS1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.