




ALDH1A2 Antibody

Product Code	CSB-PA144510
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94788
Immunogen	Synthesized peptide derived from internal of HumanALDH1A2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ALDH1A2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AL1A2; aldehyde dehydrogenase 1 family; member A2; aldehyde dehydrogenase 1A2; aldehyde dehydrogenase family 1 member A2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ALDH1A2

Image



Western blot analysis of extracts from K562 cells, using ALDH1A2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.