




ARHGEF1 Antibody

Product Code	CSB-PA791509
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92888
Immunogen	Synthesized peptide derived from internal of Human ARHGEF1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ARHGEF1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	p115-RhoGEF; p115RhoGEF; 115 kDa guanine nucleotide exchange factor; Sub1.5;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARHGEF1

Image



Western blot analysis of extracts from COLO205, HepG2, HuvEc and Jurkat cells, using ARHGEF1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.