




ARHGEF12 Antibody

Product Code	CSB-PA111426
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZN5
Immunogen	Synthesized peptide derived from internal of Human ARHGEF12.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total ARHGEF12 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ARHGC; ARHGEF12; Leukemia-associated RhoGEF; RHO guanine nucleotide exchange factor (GEF) 12 (LARG); Rho guanine nucleotide exchange factor 12
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARHGEF12

Image



Western blot analysis of extracts from RAW264.7 cells, using ARHGEF12 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.