




ARHGEF19 Antibody

Product Code	CSB-PA925991
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IW93
Immunogen	Synthesized peptide derived from internal of Human ARHGEF19.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARHGEF19

Image



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using ARHGEF19 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.