




RHOB Antibody

Product Code	CSB-PA916815
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62745
Immunogen	Synthesized peptide derived from internal of Human RHOB.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total RHOB protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	H6; Rho-related GTP-binding protein RhoB; RHOB; ARH6; ARHB
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RHOB

Image



Western blot analysis of extracts from HepG2 cells and HeLa cells, using RHOB antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.