





SEPT6 Antibody

Product Code	CSB-PA246269
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14141
Immunogen	Synthesized peptide derived from internal of Human SEPT6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SEPT6 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Septin-6; SEPT6;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SEPT6

Image



Western blot analysis of extracts from Jurkat cells, using SEPT6 antibody.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using SEPT6 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.