




SEPT7 Antibody

Product Code	CSB-PA933190
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16181
Immunogen	Synthesized peptide derived from internal of Human SEPT7.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total SEPT7 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CDC10; SEP7; septin 7;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SEPT7

Image



Western blot analysis of extracts from MCF-7 cell, HuvEc cells and HeLa cells, using SEPT7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.