




PRSS33 Antibody

Product Code	CSB-PA138994
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NF86
Immunogen	Synthesized peptide derived from internal of Human PRSS33.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PRSS33 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Serine protease 33 [Precursor]; EC 3.4.21.-; Serine protease EOS; PRSS33;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRSS33

Image



Western blot analysis of extracts from 293 cells, COLO cells and Jurkat cells, using PRSS33 antibody.