




NEK5 Antibody

Product Code	CSB-PA821453
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P3R8
Immunogen	Synthesized peptide derived from internal of Human NEK5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NEK5 protein.
Tested Applications	ELISA, WB, IHC, IF; WB:1:500-1:3000, IHC:1:50-1:200, IF:1:100-1:200
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	NimA-related protein kinase 5;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NEK5

Image



Western blot analysis of extracts from HepG2 cells and COLO205 cells, using NEK5 antibody.