




PLK3 Antibody

Product Code	CSB-PA114926
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4B4
Immunogen	Synthesized peptide derived from internal of Human PLK3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PLK3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.21; Polo-like kinase 3; PLK-3; Cytokine-inducible serine/threonine-protein kinase; FGF-inducible kinase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PLK3

Image



Western blot analysis of extracts from HepG2 cells, using PLK3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.