




PPP2R5A Antibody

Product Code	CSB-PA000119
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15172
Immunogen	Synthesized peptide derived from internal of Human PPP2R5A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PPP2R5A protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit alpha isoform; PP2A B subunit isoform B-alpha; PP2A B subunit isoform B56-alpha; PP2A B subunit isoform PR61-alpha; PR61alpha
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP2R5A

Image



Western blot analysis of extracts from HuvEc cells and HeLa cells, using PPP2R5A antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.