





PPP2R1B Antibody

Product Code	CSB-PA830420
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30154
Immunogen	Synthesized peptide derived from internal of Human PPP2R1B.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PPP2R1B protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	2AAB; beta isoform of regulatory subunit A; protein phosphatase 2; PP2A; A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP2R1B

Image



Western blot analysis of extracts from Jurkat cells and HeLa cells, using PPP2R1B antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using PPP2R1B antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.