




PPP4R1L Antibody

Product Code	CSB-PA090551
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P1A2
Immunogen	Synthesized peptide derived from internal of Human PPP4R1L.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PPP4R1L protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	F118B; FAM118B; FLJ21103;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP4R1L

Image



Western blot analysis of extracts from HepG2 cells, COLO cells, Jurkat cells, HUVEC cells and COS cells, using PPP4R1L antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.