




PPP4R2 Antibody

Product Code	CSB-PA124096
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NY27
Immunogen	Synthesized peptide derived from internal of Human PP4R2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP4R2

Image



Western blot analysis of extracts from K562 cells and COS-7 cells, using PPP4R2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.