





SERPINB9 Antibody

Product Code	CSB-PA065366
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50453
Immunogen	Synthesized peptide derived from internal of Human SERPINB9.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SERPINB9 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cytoplasmic antiproteinase 3; CAP-3; CAP3; Proteinase inhibitor 9;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SERPINB9

Image



Western blot analysis of extracts from Jurkat cells, using SERPINB9 antibody.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using SERPINB9 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.