




SHB Antibody

Product Code	CSB-PA125145
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15464
Immunogen	Synthesized peptide derived from C-terminal of Human SHB.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SHB protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	SH2 domain-containing adapter protein B; SHB;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SHB

Image



Western blot analysis of extracts from HuvEc cells, using SHB antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.