




SH3RF3 Antibody

Product Code	CSB-PA243176
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TEJ3
Immunogen	Synthesized peptide derived from C-terminal of Human SH3RF3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total SH3RF3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	FLJ00204; Plenty of SH3s 2; plenty of SH3s-2; POSH2; SH3 domain containing ring finger 3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SH3RF3

Image



Western blot analysis of extracts from HuvEc cells, using SH3RF3 antibody.