




SLC24A4 Antibody

Product Code	CSB-PA129178
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFF2
Immunogen	Synthesized peptide derived from internal of Human SLC24A4.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total SLC24A4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Na(+)/K(+)/Ca(2+)-exchange protein 4; NCKX4; sodium/potassium/calcium exchanger 4; solute carrier family 24 (sodium/potassium/calcium exchanger); mem
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC24A4

Image



Western blot analysis of extracts from RAW264.7 cells, using SLC24A4 antibody.