




ATP1B3 Antibody

Product Code	CSB-PA297832
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54709
Immunogen	Synthesized peptide derived from Internal of Human ATP1B3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ATP1B3 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AT1B3; ATP1B3; ATPase; Na ⁺ /K ⁺ transporting; beta 3 polypeptide
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ATP1B3

Image



Western blot analysis of extracts from K562 cells, using ATP1B3 antibody.