




SPTBN5 Antibody

Product Code	CSB-PA208853
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRC6
Immunogen	Synthesized peptide derived from internal of Human SPTBN5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SPTBN5 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Spectrin beta chain; brain 4; Spectrin; non-erythroid beta chain 4; Beta-V spectrin
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SPTBN5

Image



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using SPTBN5 antibody.