




SRSF11 Antibody

Product Code	CSB-PA695198
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05519
Immunogen	Synthesized peptide derived from internal of Human SFRS11.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total SFRS11 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Arginine-rich 54 kDa nuclear protein; p54; SFR11; SFRS11; Splicing factor arginine/serine-rich 11
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SRSF11

Image



Western blot analysis of extracts from LOVO cells and RAW264.7 cells, using SFRS11 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.