





TUBGCP6 Antibody

Product Code	CSB-PA905394
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RT7
Immunogen	Synthesized peptide derived from internal of Human TUBGCP6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TUBGCP6 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	gamma-tubulin complex component 6; GCP-6; GCP6; KIAA1669;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TUBGCP6

Image



Western blot analysis of extracts from COLO cells, using TUBGCP6 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using TUBGCP6 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.