




TUBGCP4 Antibody

Product Code	CSB-PA020361
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGJ1
Immunogen	Synthesized peptide derived from internal of Human TUBGCP4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total TUBGCP4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Gamma-tubulin complex component 4; GCP-4; hGCP4; Hgrip76; h76p
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TUBGCP4

Image



Western blot analysis of extracts from COLO cells and HT-29 cells, using TUBGCP4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.