




SVOP Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA056297 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8N4V2 |
| Immunogen | Synthesized peptide derived from internal of Human SVOP. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse |
| Specificity | The antibody detects endogenous levels of total SVOP protein. |
| Tested Applications | ELISA,WB;WB:1:500-1:3000 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | Synaptic vesicle 2-related protein; SV2-related protein; |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | SVOP |

Image



Western blot analysis of extracts from K562 cells, using SVOP antibody.