




SVOP Antibody

Product Code	CSB-PA056297
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N4V2
Immunogen	Synthesized peptide derived from internal of Human SVOP.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total SVOP protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Synaptic vesicle 2-related protein; SV2-related protein;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SVOP

Image



Western blot analysis of extracts from K562 cells, using SVOP antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.