




SYT16 Antibody

Product Code	CSB-PA153617
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q17RD7
Immunogen	Synthesized peptide derived from internal of Human SYT16.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SYT16 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Synaptotagmin-16; Synaptotagmin XIV-related protein; Synaptotagmin 14-like protein; Chr14Syt; SYT16
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SYT16

Image



Western blot analysis of extracts from K562 cells, Jurkat cells, 293 cells and HepG2 cells, using SYT16 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.