




SYTL4 Antibody

Product Code	CSB-PA780689
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C24
Immunogen	Synthesized peptide derived from internal of Human SYTL4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total SYTL4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	exophilin-2; granuphilin; granuphilin-a; synaptotagmin-like protein 4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SYTL4

Image



Western blot analysis of extracts from 3T3 cells, using SYTL4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.