




CD1E Antibody

Product Code	CSB-PA912791
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15812
Immunogen	Synthesized peptide derived from internal of Human CD1E.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CD1E protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	T-cell surface glycoprotein CD1e; CD1e antigen; R2G1; CD1E;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD1E

Image



Western blot analysis of extracts from HuvEc cells and HepG2 cells, using CD1E antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.