




CCT7 Antibody

Product Code	CSB-PA963179
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99832
Immunogen	Synthesized peptide derived from internal of Human CCT7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CCT7 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	TCP-1-eta; CCT-eta; HIV-1 Nef-interacting protein; CCTH; NIP7-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CCT7

Image



Western blot analysis of extracts from HeLa cells, using CCT7 antibody.