





CCT6A Antibody

Product Code	CSB-PA796607
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40227
Immunogen	Synthesized peptide derived from internal of Human CCT6A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CCT6A protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	T-complex protein 1 subunit zeta; TCP-1-zeta; CCT-zeta; CCT-zeta-1; Tcp20
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CCT6A


Image



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using CCT6A antibody.



Western blot analysis of extracts from HepG2 cells, using CCT6A antibody.



Western blot analysis of extracts from cos-7 cells (Lane 2) and 293 cells (Lane 3), using CCT6A antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.