




TSSK4 Antibody

Product Code	CSB-PA877181
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6SA08
Immunogen	Synthesized peptide derived from internal of Human TSSK4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total TSSK4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	C14orf20; EC 2.7.11.1; serine/threonine-protein kinase 22E; STK22E; testis-specific kinase 4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TSSK4

Image



Western blot analysis of extracts from HT-29 cells, using TSSK4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.