




TCEAL3/TCEAL5/TCEAL6 Antibody

Product Code	CSB-PA094518
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969E4
Immunogen	Synthesized peptide derived from internal of Human TCEAL3/5/6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TCEAL3/5/6 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transcription elongation factor A protein-like 3; TCEA-like protein 3; Transcription elongation factor S-II protein-like 3; TCEAL3; MSTP072
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TCEAL3/TCEAL5/TCEAL6

Image



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using TCEAL3/5/6 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.