




TCEAL4 Antibody

Product Code	CSB-PA136717
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EI5
Immunogen	Synthesized peptide derived from internal of Human TCEAL4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TCEAL4 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transcription elongation factor A protein-like 4; TCEA-like protein 4; Transcription elongation factor S-II protein-like 4; TCEAL4; NPD017
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TCEAL4

Image



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using TCEAL4 antibody.