




TCF7L1 Antibody

Product Code	CSB-PA737037
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCS4
Immunogen	Synthesized peptide derived from internal of Human TCF7L1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TCF7L1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transcription factor 7-like 1; HMG box transcription factor 3; TCF-3; TCF7L1; TCF3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TCF7L1

Image



Western blot analysis of extracts from 3T3 cells, using TCF7L1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.