




TAF3 Antibody

Product Code	CSB-PA564583
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VWG9
Immunogen	Synthesized peptide derived from internal of Human TAF3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total TAF3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transcription initiation factor TFIID subunit 3; 140 kDa TATA box-binding protein-associated factor; TBP-associated factor 3; Transcription initiation factor TFIID 140 kDa subunit; TAF(II)140
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TAF3

Image



Western blot analysis of extracts from K562 cells, using TAF3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.