




EIF2B3 Antibody

Product Code	CSB-PA912439
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR50
Immunogen	Synthesized peptide derived from Internal of Human EIF2B3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total EIF2B3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EI2BG; EIF-2B; eIF-2B GDP-GTP exchange factor subunit gamma; EIF2B3; EIF2Bgamma
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF2B3

Image



Western blot analysis of extracts from Jurkat cells, using EIF2B3 antibody.