





# TSEN54 Antibody

<b>Product Code</b>	CSB-PA272793
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q7Z6J9
<b>Immunogen</b>	Synthesized peptide derived from internal of Human TSEN54.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Specificity</b>	The antibody detects endogenous levels of total TSEN54 protein.
<b>Tested Applications</b>	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	tRNA-splicing endonuclease subunit Sen54; tRNA-intron endonuclease Sen54; HsSen54; TSEN54; SEN54
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	TSEN54

## Image



Western blot analysis of extracts from 3T3 cells and Jurkat cells, using TSEN54 antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using TSEN54 antibody.

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.