




PEAK1 Antibody

Product Code	CSB-PA981746
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H792
Immunogen	Synthesized peptide derived from N-terminal of Human SGK269.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total SGK269 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Pseudopodium-enriched atypical kinase 1; Sugen kinase 269; Tyrosine-protein kinase SgK269; PEAK1; KIAA2002
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PEAK1

Image



Western blot analysis of extracts from 3T3 cells, using SGK269 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.